

Testele de control pentru produsul finit:

CLORURA DE SODIU 0.9% - 250ML

Nr. Lot: 508105

Data productie: 27/08/2015

Data expirare: 27/08/2018

Test	Valoare Nominala	Conform	Metoda	Rezultat
Aspect	Solutie clara, fara culoare	pH.Eur., 2.2.1 % 2.2.2 Metoda I	Vizuala	Conform
pH	5.0-7.0	pH.Eur., 2.2.3	Potentiometrie	5.6
Sodiu	154,0 mmol/L \pm 5%	pH.Eur., 2.2.22	Flamfotometru	153.30
Clorura	154,0 mmol/L \pm 5%	USP Monografie ce contine test pt. precipitarea clorurii, prin titrare	Volumetrie	153.9
Identificare: - Sodiu - Clorura	Determinare cantitativa: Pozitiva			Conform
Sterilitate	Sterila	pH.Eur., 2.3.1	Microbiologie	Conform
Endotoxine bacteriene	<0.25 I.U/ml	pH.Eur., 2.6.1	LAL- test	Conform
Data expirare	Termen valab. 3ani			08/2018
Volum extractabil	> volumul decalarat	pH.Eur., 2.6.14		Conform
Contaminarea de particule (noninvizibile)	$\geq 5\mu\text{m} \leq 100/\text{ml}$ $\geq 10\mu\text{m} \leq 25/\text{ml}$ $\geq 25\mu\text{m} \leq 3/\text{ml}$	pH.Eur., 2.9.17	Luminozitate	3.2 0.3 0.00

Responsabil Expert:

Chimist:

SS . indescifrabil

Biolog:

S.S. indescifrabil



DATA: 15.09.2015

Responsabil controlul calitatii

Vasiliki Kefala

S.S. indescifrabil