

VIOSER S.A
QUALITY CONTROL DPT.

Control Tests on the Finished Product:

SOD. CHLORIDE 0.9% ROMANIA - 1000ml

Batch No: 704557

Prod. Date: 12/04/2017

Expiry Date: 12/04/2020

Test for	Nominal Value	According to	Method	Test Results
Appearance	clear, colorless solution	Ph.Eur.2.2.1 & 2.2.2 Method I	Visually	Conform
pH	5.0-7.0	Ph.Eur.2.2.3	potentiometer	5.67
Sodium	154.0 mmol/L ± 5%	Ph.Eur.2.2.22	flame photometer	154.20
Chloride	154.0 mmol/L ± 5%	USP Monograph containing assay for chlorides precipitation titration	Volumetric	154.80
Identification				
-Sodium	Quantitative determination serves as proof for identity			Conform
-Chloride	positive	Ph.Eur.2.3.1		Conform
Sterility	sterile	Ph.Eur.2.6.1	Microbiolog.	Conform
Bacterial Endotoxins	< 0.25 IU/ml	Ph.Eur.2.6.14	LAL-test	Conform
Expiry Date	shelf life: 3 years			04/2020
Extractable Volume	> declared volume	Ph.Eur.2.9.17		Conform
Particulate Contamination	≥ 5 µm ≤ 100/ml	Ph.Eur.2.9.19	Light Blockage	3.10
-nonvisible particles-	≥ 10 µm ≤ 25/ml			0.60
	≥ 25 µm ≤ 3/ml			0.00

EXPERT RESPONSIBLES:

DATE: 27/04/2017

Chemist: 

Quality Control Manager

Biologist: 

Vasiliki Kefala 