

**VIOSER S.A**  
**QUALITY CONTROL DPT.**

**Control Tests on the Finished Product:**

**SOD. CHLORIDE 0.9% ROMANIA - 250ml**

Batch No: 706955  
Prod. Date: 21/06/2017  
Expiry Date: 21/06/2020

Test for	Nominal Value	According to	Method	Test Results
Appearance	clear, colorless solution	Ph.Eur.2.2.1 & 2.2.2 Method I	Visually	Conform
pH	5.0-7.0	Ph.Eur.2.2.3	potentiometer	5.76
Sodium	154.0 mmol/L $\pm$ 5%	Ph.Eur.2.2.22	flame photometer	154.60
Chloride	154.0 mmol/L $\pm$ 5%	USP Monograph containing assay for chlorides precipitation titration	Volumetric	154.60
<b>Identification</b>				
-Sodium	Quantitative determination serves as proof for identity			Conform
-Chloride	positive	Ph.Eur.2.3.1		Conform
Sterility	sterile	Ph.Eur.2.6.1	Microbiolog.	Conform
Bacterial Endotoxins	< 0.25 IU/ml	Ph.Eur.2.6.14	LAL-test	Conform
Expiry Date	shelf life: 3 years			06/2020
Extractable Volume	> declared volume	Ph.Eur.2.9.17		Conform
Particulate Contamination -nonvisible particles-	$\geq 5 \mu\text{m} \leq 100/\text{ml}$ $\geq 10 \mu\text{m} \leq 25/\text{ml}$ $\geq 25 \mu\text{m} \leq 3/\text{ml}$	Ph.Eur.2.9.19	Light Blockage	0.50 0.10 0.00

**EXPERT RESPONSIBLES:**

Chemist:

Biologist:

**DATE: 07/07/2017**

**Quality Control Manager**

**Dr. Theofilos Vlassis**